


| % Undersize | 3 | 5 | 10 | 25 | 50 | 75 | 90 | 95 | 97 |
|-----------------------------|------|------|------|------|-------------|-------|-------|-------|-------|
| 1.5μm | 1.00 | 1.08 | 1.19 | 1.43 | 1.61 | 1.77 | 1.91 | 2.01 | 2.06 |
| 2μm | 1.39 | 1.49 | 1.64 | 1.86 | 2.04 | 2.22 | 2.39 | 2.51 | 2.60 |
| 3μm | 2.46 | 2.57 | 2.78 | 3.06 | 3.31 | 3.56 | 3.80 | 3.93 | 4.03 |
| 4.5μm | 4.06 | 4.15 | 4.27 | 4.45 | 4.61 | 4.75 | 4.85 | 4.94 | 4.99 |
| 5μm | 3.97 | 4.07 | 4.20 | 4.43 | 4.68 | 4.99 | 5.35 | 5.56 | 5.72 |
| 5.5μm | 4.72 | 4.84 | 5.03 | 5.28 | 5.54 | 5.87 | 6.15 | 6.30 | 6.38 |
| 7μm | 5.38 | 5.74 | 6.11 | 6.52 | 6.79 | 7.03 | 7.20 | 7.32 | 7.39 |
| 7.5μm | 6.45 | 6.56 | 6.79 | 7.17 | 7.55 | 8.08 | 8.93 | 9.59 | 9.95 |
| 9μm | 7.66 | 7.89 | 8.21 | 8.63 | 9.06 | 9.45 | 10.13 | 11.28 | 11.74 |
| 10μm | 8.21 | 8.43 | 8.74 | 9.20 | 9.64 | 10.15 | 10.91 | 11.99 | 12.78 |

(by Number) (Data provided as a guide only, microspheres un-certified)



Dr G R Rideal
Founder & Senior Analyst

MONODISPERSE SILICA PARTICLE STANDARD



Whitehouse Scientific

www.WhitehouseScientific.com

NON-CERTIFIED
MONODISPERSE SILICA STANDARD



1.5, 2, 3, 4.5, 5, 5.5, 7, 7.5, 9, 10 μ m

Catalogue No: MSS1-10 Nominal Weight: 10x 0.2g FOR REFERENCE ONLY